

タイテック製品 フロンガスリスト

製品型名	冷媒番号	備考
150L	R12	
150LF	R134a	
150LF-Titan	R134a	
150LFチタン	R134a	
150L-Titan	R12	
150W	R12	
200D	R13B1	
250DF	R507A	
300L	R12	
300LF	R22	
300LF-Titan	R22	
300L-Titan	R12	
400D	R13B1	
400H	R12	
400R	R404A	
600H	R12	
600L	R12	
600LF	R22	
600L-Titan	R12	
800D	1次側 R502 2次側 R503	
800DF	1次側 R22 2次側 R23	
80L	R12	
80LF	R134a	
80W	R12	
A- I	R22	
A- II	R12	
BC-1200	R12	
BC-2200	R12	
BC-3200	R12	
BCL-19	R12	
BCP-120F	R134a	
BCP-320F	R22	
BR-15LF	R134a	
BR-160L	R12	
BR-160LF	R134a	
BR-180LF	R134a	
BR-200L	R12	
BR-3000L	R12	
BR-3000L(光)	R12	
BR-3000LF	R134a	
BR-3000LF(光)(旧)	R134a	2002.5.21以前
BR-3000LF(光)(新)	R22	2002.5.22以後
BR-300L	R12	
BR-300L(光)	R12	
BR-300LF	R134a	
BR-300LF(光)	R134a	
BR-30L	R12	
BR-30LF	R134a	
BR-31FL	R134a	
BR-32FL	R134a	
BR-33FL	R134a	
BR-3xFL	R134a	
BR-40L	R12	
BR-40LF	R134a	
BR-4xFL	R134a	
BR-50L	R12	
CB-40	R134a	
CH-400AHS-MO-FH	R22	

CH-401AF		R22	
CH-401BF		R22	
CH-402BF		R22	
CH-601AF		R22	
CH-601BF		R22	
CH-602BF		R22	
CL-10		R12	
CL-150		R12	
CL-150F		R134a	
CL-150R		R134a	
CL-19		R12	
CL-30		R12	
CL-300		R12	
CL-300C		R12	
CL-300F		R22	
CL-300FC		R22	
CL-300R		R22	
CL-400		R12	
CL-400F		R22	
CL-600		R12	
CL-600F		R22	
CL-600R		R22	
CL-80		R12	
CL-80F		R134a	
CL-80R		R134a	
CMH-Ⅲ		R12	
CML-Ⅰ		R12	
CML-Ⅲ		R12	
CP-150		R12	
CP-150F		R134a	
CP-150R		R134a	
CP-300		R12	
CP-300F		R22	
CP-300R		R22	
CP-400		R12	
CP-400F		R22	
CP-600		R12	
CP-600F		R22	
CP-80		R12	
CP-80F		R134a	
CP-80R		R134a	
CT-400		R22	
DL-2(旧)	1次側 2次側	R502 R13	
DL-2(新)	1次側 2次側	R502 R503	
EH-1500BS		R22	
EH-1500BW		R22	
EH-2200BW(旧)		R22	2009.6.17以前
EH-3750BW(旧)		R22	2009.6.17以前
EL-15		R12	
EL-150		R12	
EL-15F		R134a	
EL-400		R12	
EL-600		R12	
EL-75		R12	
EL-8		R12	
EL-8F		R134a	
EZL-80		R12	
EZL-80F		R134a	
EZL-81F		R134a	
FC-4		R12	
FC-6		R12	

G・BR-300		R134a	
GC-80		R12	
GC-80		R12	
LB-21JR		R134a	
LFD-130		R600a	
LFD-130N		R600a	
LFD-300L		R509A	
LS- I		R12	
LS- II		R12	
LX-2100		R12	
LX-2200		R12	
LX-2300F		R134a	
LX-3000		R12	
LX-3200F		R134a	
M-100L		R12	
M-200		R12	
M-200F		R22	
M-210		R12	
M-210FN		R22	
M-230		R12	
M-230F		R134a	
M-230FN		R134a	
M-240		R12	
M-250		R12	
M-260		R502	
M-260F		R134a	
M-280		R134a	
M-510		R12	
M-600		R12	
M-600F		R22	
M-600FN(旧)		R22	2003.01まで
M-600FN(新)		R-134a	2003.02以降
ML-10		R12	
ML-10F		R134a	
N-210F		R22	
TP-130		R12	
TR-80		R-502	
VA-150		R13B1	
VA-250F		R404a	
VA-300	1次側 2次側	R502 R503	
VA-500F	1次側 2次側	R22 TP5R3	
VA-500R	1次側 2次側	R404A TP5R3	
VA-800	1次側 2次側	R12 R503	
VA-800F	1次側 2次側	R22 R23	
VA-800特	1次側 2次側	R502 R503	
VD-15		R13B1	
VD-16		R13B1	
VD-250F		R507A	
VD-30	1次側 2次側	R502 R503	
VD-31	1次側 2次側	R502 R503	
VD-400F		R22	
VD-41		R13B1	
VD-500F	1次側 2次側	R22 R508A	

VD-500R	1次側 2次側	R404A R508A	
VD-60		R13B1	
VD-600		R502	
VD-61		R13B1	
VD-800F	1次側 2次側	R22 R23	
VD-81	1次側 2次側	R12 R503	
WB-5L		R12	
YW-100		R134a	
YW-12		R134a	
YW-300		R22	
YW-75		R134a	
ZL-100		R134a	
パーソナルLt		R12	
パーソナルLt-10		R12	
パーソナルLt-10F		R134a	
メディカルサーモバス		R12	
メディカルサーモバスⅡ		R12	

2017.7.3更新